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23.11.04 CB/GP/TM

Operating Manual Mode d'emploi

Pumpen-Antrieb BVP Process (IP 65) ISM 920

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Pro 280 / Pro 281
Pro 380 / Pro 381
360 / 380 / 380 AD
CA 4 / CA 8 / CA 12
FMI Q0 / Q1 / Q2 / Q3
Easy-Load / Quickload
Standard
MS/CA 4-12 / MS/CA 8-6
PTFE Tube 2 mm / 4 mm
SB 2V / 3V
WM 5

Deutsch

Pump drive BVP Process (IP 65) ISM 920

For pump-heads

Pro 280 / Pro 281
Pro 380 / Pro 381
360 / 380 / 380 AD
CA 4 / CA 8 / CA 12
FMI Q0 / Q1 / Q2 / Q3
Easy-Load / Quickload
Standard
MS/CA 4-12 / MS/CA 8-6
PTFE Tube 2 mm / 4 mm
SB 2V / 3V
WM 5

English

Moteur BVP Process (IP 65) ISM 920

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Pro 280 / Pro 281
Pro 380 / Pro 381
360 / 380 / 380 AD
CA 4 / CA 8 / CA 12
FMI Q0 / Q1 / Q2 / Q3
Easy-Load / Quickload
Standard
MS/CA 4-12 / MS/CA 8-6
PTFE Tube 2 mm / 4 mm
SB 2V / 3V
WM 5

Français

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Contents

Safety precautions	4–7
Warranty terms	8–9
Product	10
Rear panel	11
Mains voltage	12
Starting the pump	13
Operating panel	14
Start-up information	15
Control keys	16
Pumping	17
Pumping against pressure	18
When the pump is not in use	19
Tubing	20
Running-in period	20
Tubing life	20
Analog interface	21–23
Switch S1 and S2	24
Access to the switches	25

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Contents

Accessories	
Foot switch	26
Spare-cassettes	26–27
Interchangeable Pump-heads	28–29
Flow rates charts	30–35
Opening the casing hood	36–38
Plug gaskets	39
Changing the fuses	40
Maintenance	41
Repairs	41
Spare-parts	42
Technical specifications	42–43

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Please note

We recommend you to read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

ISMATEC SA does not take liability for any damage resulting from the use of an ISMATEC® pump.

ISMATEC SA does not admit responsibility for the handling of chemicals.

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Safety precautions

Ismatec® Pumps are designed for pumping applications in laboratories and industry.

As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed.

- The circuit between mains supply and pump has to be earthed.
- The pump must not be operated outside the destined operating and environmental conditions.
- The pump must not be used:
 - for medical applications on human beings
 - in explosion proof chambers or in the presence of flammable gases or fumes

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Safety precautions

- The pump must be switched off when pump-heads, cassettes or tubing are inserted or changed.
- The permeability of tubing depends on the material used and pressure conditions. Tubing can also become electro-statically charged. Please be aware of possible hazards when laying tubing in explosion-proof chambers.
- Pump-heads consist of revolving parts. Therefore, the pump must not be operated before the cassettes are fully snapped-in or the tube-bed completely shut.
- Do not manipulate the pump-head before the pump is switched off and disconnected from the mains supply.

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Safety precautions

- Be particularly cautious that no parts of your body such as fingers, long hair, etc. or jewellery, or loose objects such as cables or tubing, etc. can be trapped by the revolving rotor.
- Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken. When pumping corrosive media we recommend to place the pump in a collecting basin.
- Only new fuses, according to the specifications stated on page 40 in this manual, must be used.
- Do not open or remove the housing while the pump is operating.

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Safety precautions

- Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- For service and repairs carried out by the customer or by third-party companies Ismatec SA denies any responsibility.

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✉ Warranty

For all parts manufactured by Ismatec® from date of delivery: 2 years

All other parts, excluding consumables: 1 year

In case of any queries, please contact your Ismatec® representative.

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Warranty terms

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion.

A defective pump must be returned in the original Ismatec® packing or in a packet of equal quality. For pump-heads from manufacturers other than Ismatec SA the warranty terms of the specific manufacturer are valid.

The duration of the warranty is not affected by making a claim for warranty service.

Further claims are excluded. Shipping costs are charged to the customer.

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Warranty terms

Our warranty becomes invalid in the case of:

- improper operation by the user, or if the pump is diverted from its proper use
- unauthorized modification or misuse by the user or by a third party
- improper site preparation and maintenance
- operation outside of the environmental and electrical specifications for the product
- use of third-party software, hardware, accessories or consumables purchased by the user and which do not comply with our specifications.

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Product

Contents of the package

- BVP *Process* drive
Order No. ISM 920
incl. integrated power cord,
2 m long, with IEC 320 plug (male)
- 1 power cord with an IEC 320 con-
nector (female plug) and a country
specific mains plug
- 1 gaskets (in plastic bag) for
analog interface (AD0069)
- Operating manual

If ordered:

- pump-head with mounting plate
- other accessories

Please check the package

and its contents for transport dama-
ge. If you find any signs of damage,
please contact your local Ismatec®
representative immediately.
Complaints can only be accepted
within 8 days from receipt of the
goods.

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Rear panel

- 1 Mains switch (on/off)
- 2 Power cord
- 3 Analog interface
 - input for:
 - speed control
0–5 V or 0–10 V, and
0–20 mA or 4–20mA
 - rotation direction
 - Run/Stop
 - Foot switch

- Output for:
 - speed 0–10 V_{DC}
or 0–7.2 kHz


⚠ For using the analog interface under IP-65 conditions please refer to page 39.

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* See also page 40

Socket/Power cord

The socket must be earthed. (protec-
tive conductor contact)

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Starting the pump

- Mount the pump-head according to the mounting instruction manual supplied with the pump-head
- If using an FMI piston pump-head, adjust the angle for the piston stroke volume (for detailed information please refer to the instruction manual of the pump-head)
- Insert the tubing into the peristaltic pump-head
- Connect the pump tubing to the system
- Connect the pump to the mains and switch it on with the power supply switch

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Operating panel

- 1 Digital LED display for rotation speed 1.0–240.0 rpm, adjustable in steps of 0.1 rpm
- 2 Reduce speed
- 3 Increase speed
- 4 Run / Stop
- 5 MAX key for max. speed (ideal for filling or emptying the tubing)
- 6 MIN key for min. speed (for mounting the pump-head)
- 7 Rotation direction, clockwise
- 8 Rotation direction, counter-clockwise

9 LED indicators

- Pump = pump is controlled via operating panel (blinks when pump is operating)
- Remote = pump is controlled via analog interface

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Start-up information

After switching on the power supply switch, the following values are displayed:

- 1 LED test »8.8.8.8.«
- 2 Firmware version

Only for FMI pump-heads

In order to realize the possible flow rates, the required angle must be set on the pump-head.

We advise you to check the table on page 35 and to follow the instructions stated in the separate manual for the pump-head.

Do not pump media containing particles that exceed a diameter of 0.8 mm.

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max. speed
(when pump is running)

d MIN

min. speed (when pump is idle)
turns slowly, e.g. for moun-
ting a pump-head)

e Rotation speed

- △ = increase speed
 - ▽ = decrease speed
- Maintaining pressure on
△ or ▽ key changes display
read-out into fast mode.

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Pumping

Speed 1.0–240.0 rpm, adjustable in steps of 0.1 rpm

- Enter required speed with Δ/∇ keys
• Start pump with RUN/STOP
→ The speed can also be adjusted while the pump is running.
- Maintain pressure on MAX/MIN-key while the pump is running
→ Pump is running at maximum speed (ideal for fast filling or emptying the system)
- Maintain pressure on MAX/MIN-key when pump is idle (power switched on!)
→ pump is turning at 10 rpm (for mounting the pump-head)

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Pumping against pressure

The BVP *Process* can be used for continuous duty at a differential pressure of max. 1.5 bar. For short-time operation higher differential pressures may be managed depending on the mounted pump-head and tubing i.d. (smaller i.d. is preferable).

Depending on the application conditions the FMI pump-heads allow operations at differential pressures of up to 7 bar.


In case of any doubts please contact your Ismatec® agent.

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
When the pump is not in use

When the pump is idle, we recommend you to release the tubing from pressure.

This helps to protect the tubing from unnecessary strain and prolongs its service-life.

⚠ Syphoning effect

When the tubing is released from squeezing, the fluid can flow back to the reservoir.

- a Tubing used with cassettes do not need to be removed. Releasing the cassette on one side is sufficient.
- b The single-channel pump-heads 360, 380 and Easy Load allow opening the tube-bed in order to release the tubing pressure.

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Please note

Please refer to our detailed tubing documentation.

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Running-in period for tubing

Every new tube requires a running-in period. If constant and reproducible flow rates are required, we recommend you to run new tubing in with water or the medium to be pumped for at least 1 to 3 hours before you start the application.

Tubing life

The service-life of the tubing depends on the application and the tubing material used.

Example: chemicals, rotation speed, differential pressure, temperature, viscosity, pressure on tubing, etc.

General information on the service-life is stated in our Ismatec® catalog (without obligation!).

In order to improve the lubrication and service-life of the tubing, we recommend users to spray both the tubing and the pump rollers with our silicone oil spray (Order No. SC0179).

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
Digitale Eingänge (TTL-Pegel)

Digital inputs (TTL-level)

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Pin 2, remote / Pin 3, start
Pin 4, direction / Pin 13, speed intern

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Analog interface

Pin 1, GND (ground)

Reference potential for all other inputs

Pin 2, remote

For changing between manual control and analog interface. For activating the analog interface, pin 2 must be connected with pin 1 (GND).

Pin 3, start

- In normal operation (pin 2 open) pin 3 operates the foot-switch (see DIP switch settings, page 24–25)
- In remote operation (pin 2 to GND) the pump starts when connected to pin 1 (GND)

Pin 4, direction

In the open position the pump turns clockwise; when connected to pin 1 (GND) it turns counter-clockwise.

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
Analog interface

Pin 5, speed IN

For external speed control (0–5V, 0–10V, 0–20mA, 4–20mA)
Alternatives by means of DIP switches inside the pump (see page 24–25)

Pin 7, +26V_{DC}

About +26 V_{DC} are available (max. current 1A).

Pin 9, speed OUT

The default setting is 0–10 V_{DC}, proportionally to the speed 0–240.0 rpm.
Alternatively a frequency range from 0–7.2 kHz is available. Alternatives with sliding switch S2 inside the pump (see page 24–25)

Pin 10, +5 V_{DC}

About +5 V_{DC} are available. (max. current 0.5 A)

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Analog interface

Pin 13, speed intern

Depending on how the pump is operated, pin 13 has different functions:

Analog interface not activated (Normal operation, i.e. pin 2 is open) Pin 13 serves as **auto-start function**. If pin 13 is connected to pin 1 (GND), the pump can be started and stopped directly from the power supply (the power switch must be on).

Analog interface activated (Pin 2 on GND)

- Pin 13 open:
The rotation speed is adjusted via pin 5 (speed IN).
- Pin 13 on GND:
The rotation speed can be adjusted by the speed selector on the control panel of the pump.

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
Settings of switch S1

Operation via foot-switch

With DIP switch 5 the user can switch between two possibilities:

- »FS toggle« (On/Off)
- »FS direct«
(On = as long as pressure is maintained on the foot-switch)

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Pins		DIP-Switch 1	DIP-Switch 2	DIP-Switch 3	DIP-Switch 4	DIP-Switch 5	DIP-Switch 6
Pin 3 (Foot switch)	FS toggle					OFF*	
	FS direct					ON	
Pin 5 speed IN	0–5V		OFF*	OFF*	OFF*		
	0–10V		OFF	OFF	ON		
	0–20mA		OFF	ON	OFF		
	4–20mA		ON	ON	OFF		

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Access to the switches

In order to maintain the optimal casing seal, the pump should not be opened without a particular reason. This is, however, necessary in order to gain access to the switches.

For this purpose we recommend that you proceed exactly as indicated on pages 36–39.

→ **Make sure that the pump is disconnected from the mains supply.**

⚠ The instrument should only be opened by a qualified technician! Capacitors inside the pump may still be charged even though the mains plug has been disconnected some time ago.

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Foot switch

Order No. IS 10039

This foot switch serves as a start/stop device. It is very useful when using the pump as a dispenser for filling tubes, bottles, etc.. Both hands are free for handling the bottles and tubing.

Spare-cassettes in POM

MS/CA Click'n'go

Order No. IS 3510

CA Click'n'go

Order No. IS 3710

MS/CA pressure lever*

Order No. IS 0649

CA pressure lever*

Order No. IS 0122

Adaptor for type CA IS 0123
(order 2 adaptors per cassette)

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MS/CA Click'n'go Cassette
MS/CA Click'n'go cassette

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MS/CA pressure lever*
Order No. IS 3629

CA pressure lever*
Order No. IS 3820

Adaptor for type CA IS 3861
(order 2 adaptors per cassette)

* The cassettes with pressure lever are available on request. Depending on the application, this type of cassette may provide better results at elevated differential pressure conditions.

→ Click'n go Cassettes

(The pump does not deliver the liquid)
When using new tubing for the first time, it may occur that, depending on the tubing used (hardness and diameter), the pump cannot be primed and, hence, does not deliver the liquid. If that is the case we recommend you to prime the tubing and to run the pump with the tubing inserted for about 15 to 30 minutes.

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Mounting

Each pump-head is supplied with an instruction manual.

For selecting the correct tubing, please refer to our detailed tubing documentation.

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

0.44–2800 ml/min
380
ISM 718

Interchangeable Pump-heads

These interchangeable pump-heads can also be ordered separately and mounted on the BVP *Process* drive.

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
0.002–230 ml/min
CA 4 / ISM 721
CA 8 / ISM 732
CA 12 / ISM 733

Interchangeable Pump-heads

→ for the **piston pump heads** the suffix „-W“ resp. „-WT“ means: with isolation (W) or temperature gland for crystallizing media.

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Flow rates

⚠ In the tables listed subsequently the values indicated are only approximate and determined as follows:

ml/min, per channel, with water and Tygon® tubing, without differential pressure

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⚠ Please note

We will be pleased to help the user to select the tubing material. However, the user himself has the final responsibility for the selection of the correct tubing material.

For reproducible results we recommend you to use tubing from Ismatec®

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⚠ Piston pump-heads

After use, the piston pump-heads require thorough flushing

- especially after pumping media containing salt, protein or particles.

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Opening the casing hood

⚠ In order to maintain the best possible seal of the casing, the pump should only be opened for replacing the fuses or changing the DIP-switches. We recommend to proceed strictly as follows:

Removing the casing hood

1. Remove the 4 screws with a Phillips screw driver (size 2) on each side of the casing hood.
2. Spread **slightly** the casing hood at the bottom and remove it by carefully pulling it upwards.

Mounting the casing hood

1. Before replacing the casing hood we recommend you to lubricate the sealing on the contact surface (e.g. with silicone oil spray - see page 20). In case that the two sealing ends have come apart, they must be stuck together again with Silicone sealing com-

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Opening the casing hood

2. Spread the casing hood slightly at the bottom and push it carefully over the pump. Be careful that no cables are trapped.
3. Press the casing hood downwards and fasten the screws evenly from top to bottom. When tightening the screws, be careful that the profile sealing is not trapped. If the sealing should be damaged, it must be replaced, otherwise the seal may no longer be tight.

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Opening the casing hood

Sealing the casing hood

When the casing hood fits tight, each single screw must be removed again and stopped with silicone sealing compound.

Then, the screws must be refastened again evenly from top to bottom.

→ Remove the screw - fill the hole with silicone - fasten the screw again - remove the last bit of silicone with a rag. If necessary clean up with acetone.

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

Plug gaskets

If the interface port is used under **IP-65 conditions**, the originally supplied gasket must be inserted into the male plug of the corresponding connecting cables.

- 1 Gasket for plug (15 pin, male) for analog interface (female) Spare part No. AD0069


⚠ Choosing the correct interface cable (sealed plug, e.g. moulded plug) according to the IP-rating used, is entirely the user's responsibility.

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
Changing the fuses

The fuses are fixed to the control board (above) and the power supply print (below) as illustrated opposite.

Fuses control board

- a 4.0 A, slow-blow
- b 1.6 A, fast-blow

Fuses power supply

- c 3.15 A, slow-blow

⚠ Switch the pump off, pull out the mains plug.

⚠ Remove the casing hood.
For removing and remounting the casing hood, please refer to pages 36–39.

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Service and repairs

Provided the *BVP Process* is operated properly and in compliance with this manual, the tubing is the only part that is subject to wear and tear.

For repairs send the complete pump with detailed description of the failure to your ISMATEC® distributor. Please use the original ISMATEC® packing or a packet of equal quality.

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Replacement parts

For repairs you intend to carry out by yourself (out of the warranty period), you can ask your ISMATEC® distributor for:

- parts lists
- replacement parts for pump drive
- wiring diagrams

Please give information on defect, date of purchase, serial-no., and model.

Technical Specifications

Drive

Motor type DC motor
Speed range 1.0–240.0 rpm
digitally adjustable
in steps of 0.1 rpm

Differential pressure

1.0–1.5 bar (depending on pump head) for short-time use and with small tubing i.d.'s higher pressures are possible.

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△ Please note

We also recommend you to observe our Warranty Terms as well as our Terms and Conditions of Sale.

In case of any queries, please contact your local Ismatec® representative.

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IP 65

Operating conditions

Temperature +5 to +40°C,
Rel. humidity max. 80%
not condensing, at normal environmental conditions

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
ISMATEC® - Your competent partner for demanding metering and dispensing applications.

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A Unit of IDEX Corporation
Laboratoriumstechnik GmbH