ISMATEC[®]
A Unit of IDEX Corporation

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CE 20.12.04 CB/GP/TM

Operating Manual Mode d'emploi

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BVP-Z
ISM 446

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Deutsch

Gear pump drive

with variable speed drive

BVP-Z
ISM 446

Drive without pump-head

English

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BVP-Z
ISM 446

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Français

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Contents

Safety precautions	3–5
Warranty terms	6–7
Product	8
Operating panel	9
Rear panel	9
Mains voltage	10–11
Voltage setting	11
Changing the fuses	11
Starting the pump	12–13
Analog interface	14–16
Pumping against pressure	17
Pump-head with bypass	17
Before switching the pump off	17
Pump-head	18
Mounting the pump-head	19–20
Information on pump-heads	21
Accessories	22–24
Maintenance, repairs	25
Spare-parts	26
Technical specifications	26

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Please note

We recommend you to read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

ISMATEC SA does not take liability for any damage resulting from the use of an ISMATEC® pump.

ISMATEC SA does not admit responsibility for the handling of chemicals.

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Safety precautions

ISMATEC® gear pumps are designed for pumping applications in laboratories and industry. As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed.

- The circuit between mains supply and pump has to be earthed.
- The pump must not be operated outside the destined operating and environmental conditions.
- Place the unit in a well ventilated position at least 10 cm away from walls, partitions etc. Ensure that curtains and similar materials do not cover the ventilation slits.
- The pump **must not** be used:
 - for medical applications on human beings
 - in explosion-proof chambers or in the presence of flammable gases or fumes.

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Safety precautions

- The pump must be switched off when changing the pump-head.
- The permeability of pump-tubing depends on the material used and pressure conditions. Tubing can also become electro-statically charged. Please be aware of possible hazards when laying tubing in explosion-proof chambers.
- Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken.
- Do not manipulate the pump-head before the pump is switched off and disconnected from the mains.
- Before initial start-up, the pump-head must be filled with the pumping liquid or water. Never run the pump-head under dry conditions.

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Please note

Never pump media containing particles.

If necessary, filter the media before starting your application.

Never run the pump-head under dry conditions.

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Safety precautions

- We advise you to never pump liquid containing particles.
Even the smallest particles have an abrasive effect, which can ruin the gears in the course of a few seconds.
- Only new fuses, according to the specifications stated on page 10 in this manual, must be used.
- The fuse-holder must not be short-circuited.
- Do not open or remove the housing while the pump is operating.
- Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- For service and repairs carried out by the customer or by third-party companies ISMATEC SA denies any responsibility.

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From date of delivery

BVP-Z drive: 2 years
Pump-head: 1 year

In case of any queries, please contact your nearest ISMATEC® representative.

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Warranty Terms

We warrant the perfect functioning of the BVP-Z drive, provided it has been installed and operated correctly according to our operating instructions.

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion.

The pump-heads are covered by the specific terms of Micropump Ltd. Gears and seals are regarded as consumables and as such they are not covered by the warranty.

The duration of the warranty is not affected by making a claim for warranty service.

Further claims are excluded. Shipping costs for instruments under warranty are charged to the customer.

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Warranty Terms

Our warranty becomes invalid in the case of

- improper operation by the user, or if the pump is diverted from its proper use,
- unauthorized modification or misuse by the user or by a third party,
- improper site preparation and maintenance,
- operation outside environmental and electrical specifications for the product,
- use of third-party software, hardware, interfacing or consumables purchased by the user and which do not comply with our specifications,
- damages caused by contamination or leaks due to torn or bursted tubing.

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Product

Package contents

- BVP-Z drive
- Power cord
- Tubing clamps (2 pcs.)
- Non-return valve (Polypropylene)
- Operating manual

Accessories as ordered

- Pump head with
1 magnet, 2 tubing nipples

Please check the package and its contents for transport damage. If you find any signs of damage, please contact your local ISIMATEC® representative immediately.

Complaints can only be accepted within 8 days from receipt of the goods.

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
Operating panel

- 1 Mains switch
- 2 Control light
- 3 Speed selector
999 = 6000 rpm (max. speed)
= 100 rps
adjustable in 0.1% steps
- 4 RUN/STOP
- 5 MAX key for drive speed
(ideal for fast filling or emptying the
system)

⚠ The drive runs only in
clock-wise direction.

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

*slow-blow

⚠ Before starting-up

Check if the voltage setting visible in the window of the fuse-holder complies with your local mains voltage.

If necessary, the voltage setting must be changed and the 2 fuses must be replaced.

⚠ Socket/Power cord

Use exclusively the original power cord supplied with the pump. The socket must be earthed (protective conductor contact).

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
⚡ Mains voltage

Voltage setting 115V / 230V and changing the fuses

➔ Switch the pump off,
pull out the mains plug.

- 1 Pull out the fuse-holder by opening it with a small screw-driver (size 0).
- 2 Take out the voltage selector plate. Turn it and re-insert it into the fuse-holder so that the required voltage rating is facing the window of the fuse-holder.
- 3 Insert 2 new fuses
➔ 230 V: 2 x 1.25 A (slow-blow)
➔ 115 V: 2 x 2.50 A (slow-blow)
⚠ Use always 2 slow-blow fuses of the same type complying with the local mains voltage. In case of a defective fuse, always replace both.
- 4 Shut the fuse-holder. The voltage rating is visible in the window.

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⚠ Please note

Before initial start-up, the pump-head must be filled with the medium or water. Never run the pump-head under dry conditions. Never pump liquids that contain particles. If necessary, filter the medium before use.

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
Starting the pump

- Check the voltage rating in the window of the fuse-holder (on rear panel). If necessary, change the rating as indicated on page 11.
- If necessary, mount the pump-head as stated on pages 19–20.
- Before screwing the tubing nipples into the pump-head, we recommend you to wrap PTFE tape around the thread of the nipple. This ensures optimum sealing and helps to protect the thread.
- Mount the tubing and connect the hose-line to your system. Fasten the tubing to the nipples by means of the tubing clamps.
- Connect the pump to the mains.
- Fill the pump-head with liquid before you start pumping.

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Starting the pump

- 1 Press in the mains switch
- 2 The control light is on
- 3 Set the required speed on the 3-digit speed selector.
(999 = 6000 rpm)
→ The rotation speed can be adjusted while the pump is running.
- 4 Start the pump by switching the RUN/STOP button.
- 5 The MAX- button allows rapid filling (or emptying) of the system while the pump is running.

⚠ Before you start pumping, fill the pump-head with the medium or water. Never run the pump-head under dry conditions.

⚠ Never pump liquids that contain particles. If necessary, filter the medium before use.

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Overcurrent protector

The tubing pump BVP *Standard* features an overload protector.

When the overcurrent protector is activated, the pump is stopped and the control light flashes (see also the description of pin 12 of analog interface on page 19).

Whenever this situation occurs, **the pump must be switched off immediately.**

→ Let the pump cool down (it takes about 2 minutes until the pump is ready again)

Before the pump is re-started, it is most important to check the reason for the overload (eg. too high differential pressure).

→ Only when the cause of the overload has been detected and the failure corrected accordingly may the pump be started again.

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
Analog interface

Pin 1, (GND) Ground

Reference potential for all other inputs.

Pin 2, remote

For changing between manual control and analog interface. For activating the analog interface, pin 2 must be connected with pin 1 (GND).

Pin 3, start

In remote operation (pin 2 to GND) the pump starts when connected to pin 1 (GND)

Pin 5, speed IN

External speed control (0–5V, 0–10V, 0–20mA, 4–20mA). Input impedance and input range can be selected via Dip-switch inside the pump (see page 16).

Pin 7, +36 V_{DC}

About +36 V_{DC} (unstable) are available the pump motor (max. load 1 A).

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
Analog interface

Pin 9, speed OUT

The default setting is 0–10 V_{DC}, proportional to the motor speed 0–6000 rpm. Alternatively a frequency range from 0 to 12 kHz is available. The input range can be selected via a dip-switch inside the pump (see page 16).

Pin 10, +5V_{DC}

This pin is connected to the stabilised +5V_{DC} voltage supply of the pump.

Pin 12, idle, overload

This pin is at +5V_{DC} while the pump runs. When idle or in overload, the voltage changes to 0V_{DC}. → Switch off the mains voltage and let the device cool down for 2 minutes.

Pin 13, internal speed

If pin 13 is connected with pin 1 (ground), the pump speed can be adjusted via the speed selector on the operating panel instead of using the signal at pin 5.

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
Pin 9 Ausgang/Output/Sortie
 A = 0–10 V_{DC} (Standard)
 B = 0–12 kHz

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Pins		DIP-Switch 1	DIP-Switch 2	DIP-Switch 3	DIP-Switch 4	DIP-Switch 5	DIP-Switch 6
Pin 5 speed IN	0–5V	20k	ON*	OFF*	OFF*	OFF*	ON*
	0–10V	38k	OFF	ON	OFF	OFF	ON
	0–20mA	130	OFF	OFF	ON	OFF	ON
	4–20mA	130	OFF	OFF	OFF	ON	OFF

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→ **Make sure that the pump is disconnected from the mains supply.**

⚠ The instrument should only be opened by a qualified technician!
 Capacitors inside the pump may still be charged even though the mains plug has been disconnected some time ago.

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Pumping against pressure

The max. differential pressure depends on the mounted gear pump-head model (see page 21).

Pump-head with bypass

Some gear pump-head models have an internal fluid bypass. With the adjusting screw the required differential pressure can be limited to a value between 0.7 bar (10 psi) and the maximum permitted differential pressure of the pump-head used.

Before switching the pump off

**⚠ If no non-return valve is moun-
ted on the discharge side of the
pump, a siphoning may occur at the
moment the pump drive is stopped!**

If the pump will not be used for some time, it is strongly recommended to thoroughly rinse the pump-head with water.

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Pump-head

In general, the BVP-Z gear pump is supplied as a complete unit, consisting of the BVP-Z drive, and as ordered:

- 1 pump-head
- 2 driving magnet
- 3 tubing nipples (2 pcs.)

For the BVP-Z, a wide range of different gear pump-heads with PTFE, PPS, PEEK, or graphite gears are available.

Refer to the table on page 21 in order to obtain more detailed information on the pump-heads available.

⚠ Before mounting a new pump-head, make sure that the correct magnet is fitted to the drive.

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Mounting the pump-head

Fig. 1

The following tools are required for mounting a pump-head:

- Allen key $3/32''$
- Phillips screwdriver (Size 0)

Remove all screws from the fixing case (a) of the drive.
(Together with the drive, two long screws are additionally supplied for fastening the large pump-heads (at the top and the bottom of the pump-head face).

Fig. 2

Fasten the driving magnet on the motor shaft (through the hole underneath).

- ➔ Ensure that there is a gap of about 3 mm between the magnet and the brim of the fixing case.

⚠ During operation the back of the magnet must not touch the housing of the pump drive.

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Fig. 1

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
Mounting the pump-head

Fig. 3

Screw the two nipples into the pump-head.

IMPORTANT

Before screwing the tubing nipples into the pump-head, we recommend **you to wrap PTFE tape** around the thread of the nipple. This ensures optimum sealing and helps to protect the thread.

Fig. 4

Screw the pump-head to the fixing-case of the drive. The face of the pump-head with the engraved model number must be directed to the user. The engraved arrow indicates the flow direction.

⚠ Fill the pump-head with liquid before the initial start-up. Avoid dry-running of the pump-head.

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Fig. 3

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For temperatures above 54°C (130°F) consult your ISMATEC distributor for proper seal selection.

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The bypass should only be used for safety reasons and not for pressure control.

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Accessories

Tubing Nipple 01 for all gear pump-heads

1/8" thread NPT
for tubing i.d. 6–7 mm (1/4")
(always order 2 pcs.)

Order-No. AR 0001

Foot switch

Order No. ISM 891

This foot switch serves as a start/stop device. It is very useful when using the pump as a dispenser for filling tubes, bottles, etc.. Both hands are free for handling the bottles and tubing.

Non-return valve (Polypropylene) (included with the BVP-Z-Drive)

Order No. SD 0007

Tubing clamps

included with the BVP-Z-Drive (2 pcs.)
Order No. BN 1001

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Maintenance

We strongly recommend to thoroughly rinse the pump-head with water after use.

Provided the BVP-Z drive is operated properly and in compliance with this manual, the moving parts such as gears, etc. are the only parts that are subject to wear and tear.

The Service Kit corresponding to the individual pump-head allows the user to replace the parts that are subject to wear and tear by himself (see page 23).

Repairs

For repairs please send the defective BVP-Z drive to your ISMATEC® representative.

Please give information on defect, date of purchase, serial-no., model.

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Spare-parts

After the warranty period your ISMATEC® representative will be pleased to send you on request:

- spare-parts
- parts lists
- wiring diagrams

Please give information on defect, date of purchase, serial-no., and model.

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△ Please note

We also recommend you to observe our Warranty Terms as well as our Terms and Conditions of Sale.

In case of any queries, please contact your local ISMATEC® representative.

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Temperature +5 to +40°C,
(+41 to 104°F)

Rel. humidity max. 80%

not condensing, at normal laboratory conditions

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
Ask for the latest
ISMATEC-catalogue.

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ISMATEC® - Your competent partner for demanding metering and dispensing applications.

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ISMATEC Laboratoriumstechnik GmbH