MCP-Z Standard

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**Gear
pump
drive**

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ISM 405

ISM 405

ISM 405

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Microprocessor
controlled

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Schutzgrad
IP 30

Protection rating
IP 30

Classe de protection
IP 30

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Drive without pump-
head

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Deutsch

English

Français

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⚠ Please note

We recommend you to read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

ISMATEC SA does not take liability for any damage resulting from the use of an ISMATEC® pump.

ISMATEC SA does not admit responsibility for the handling of chemicals.

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Safety precautions

ISMATEC® gear pumps are designed for pumping applications in laboratories and industry. As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed.

- The circuit between mains supply and pump has to be earthed.
- The pump must not be operated outside the destined operating and environmental conditions.
- Place the unit in a well ventilated position at least 10 cm away from walls, partitions etc. Ensure that curtains and similar materials do not cover the ventilation slits.
- The pump must not be used:
 - for medical applications on human beings
 - in explosion proof chambers or in the presence of flammable gases or fumes
- The pump must be switched off when changing the pump-head.
- The permeability of tubing depends on the material used and pressure conditions. Tubing can also become electro-statically charged. Please be aware of possible hazards when laying tubing in explosion-proof chambers.

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Please note

Never pump media containing particles.

If necessary, filter the media before starting your application.

Never run the pump-head under dry conditions.

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Safety precautions

- Do not manipulate the pump-head before the pump is switched off and disconnected from the mains.
- Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken.
- We advise you to never pump liquid containing particles. Even the smallest particles have an abrasive effect, which can ruin the gears in the course of a few seconds.
- Before initial start-up, the pump-head must be filled with the pumping liquid or water. Never run the pump-head under dry conditions.
- Only new fuses, according to the specifications stated on page 9–10 in this manual, must be used.
- The fuse-holder must not be short-circuited.
- Do not open or remove the housing while the pump is operating.
- Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- For service and repairs carried out by the customer or by third-party companies ISMATEC SA denies any responsibility.

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Pump drive: 2 years

Pump-head: 1 year

The pump-heads are covered by the specific terms of Micropump Ltd. Gears and seals are regarded as consumables and as such they are not covered by the warranty.

In case of any queries, please contact your nearest ISMATEC® representative.

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Warranty terms

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion.

A defective pump must be returned in the original ISMATEC® packing or in a packet of equal quality.

The duration of the warranty is not affected by making a claim for warranty service.

Further claims are excluded.

Shipping costs are charged to the customer.

Our warranty becomes invalid in the case of:

- improper operation by the user, or if the pump is diverted from its proper use
- unauthorized modification or misuse by the user or by a third party
- improper site preparation and maintenance
- operation outside of the environmental and electrical specifications for the product
- use of third-party software, hardware, accessories or consumables purchased by the user and which do not comply with our specifications.

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Product

Contents of the package

- MCP-Z *Standard* drive
Order No. ISM 405
- 1 power cord with an IEC 320
connector (female plug) and a
country specific mains plug
- 2 tubing clamps
- 1 non-return valve (PP)
- Operating manual

If ordered:

- Pump-head with
1 magnet, 2 tubing adapters

Please check the package

and its contents for transport
damage. If you find any signs of
damage, please contact your
local ISMATEC® representative
immediately.

➔ Complaints can only be
accepted within 8 days from
receipt of the goods.

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
Rear panel

- 1 Socket for foot switch and hand dispenser
- 2 Ventilator
- 3 Fuse-holder with voltage selector 115/230 V
- 4 Mains socket
- 5 RS232 In (female)
- 6 RS232 Out (male)
- 7 Valve connector
- 8 Analog interface
 - input for:
 - speed control
 - 0–5 V or 0–10 V, and 0–20 mA or 4–20mA
 - Run/Stop
 - Output for:
 - speed 0–10 V_{DC} or 0–12 kHz

Mains voltage

Mains voltage	Voltage setting	Fuse rating (slow-blow)
220–240 VAC	230 V 50/60 Hz	2 x 1.25 A
110–120 VAC	115 V 50/60 Hz	2 x 2.50 A

⚠ Before starting-up

Check if the voltage setting visible in the window of the fuse-holder complies with your local mains voltage.

If necessary, the voltage setting must be changed and the 2 fuses must be replaced.

⚠ Socket/Power cord

Use exclusively the original power cord supplied with the pump. The socket must be earthed (protective conductor contact).

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

Voltage setting 115/230 V Changing the fuses

- ➔ Switch the pump off, pull out the mains plug.
- 1 Pull out the fuse-holder by opening it with a small screwdriver (size 0).
 - 2 Take out the voltage selector plate. Turn it and re-insert it into the fuse-holder so that the required voltage rating is facing the window of the fuse-holder.
 - 3 Insert 2 new fuses
➔ 230 V: 2x1.25 A (slow-blow)
➔ 115 V: 2x2.50 A (slow-blow)

⚠ Use always 2 slow-blow fuses of the same type complying with the local mains voltage.

- 4 Shut the fuse-holder. The voltage rating is visible in the window.

Changing the fuses on the control board

The 2 fuses are fixed to the control board as illustrated opposite.

- a 5.0 A, slow-blow
- b 1.6 A, fast-blow

- ➔ Make sure that the pump is disconnected from the mains supply.
⚠ **The instrument should only be opened by a qualified technician!** Capacitors inside the pump may still be charged even though the mains plug has been disconnected some time ago.

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Starting the pump

- Check the voltage rating in the window of the fuse-holder (on rear panel). If necessary, change the rating as indicated on page 9.
- Mount the pump-head as stated on pages 42–44.
- Before screwing the tubing adapters into the pump-head, we recommend you to wrap **PTFE tape around the thread of the adapter**. This ensures optimum sealing and helps to protect the thread.
- Connect the tubing to your system. Fasten the tubing to the adapters by means of the tubing clamps.
- Connect the pump to the mains and switch it on with the power supply switch.
- ➔ Enter the ID Codes of the mounted pump-heads in the program memory currently used (see pages 15/44).

⚠ **Before you start pumping, fill the pump-head with the medium or water. Never run the pump-head under dry conditions.**

⚠ **Never pump liquids that contain particles. If necessary, filter the medium before use.**

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
Bedienungspanel

PUMP	<input type="radio"/> rpm	<input type="radio"/> Flow rate
DISP	<input type="radio"/> Time	<input type="radio"/> Volume
PAUSE	<input type="radio"/> Time	<input type="radio"/> PROGRAM
TOTAL	<input type="radio"/>	

Anzeige der Betriebsart
LEDs for operating modes
Affichage du mode
d'exploitation

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Operating panel

- 1 Mains switch (on/off)
- 2 Digital LED display
- 3 Reduce value
- 4 Increase value
- 5 Run / Stop
- 6 Calibration / Data saving
- 7 MAX/RESET
(Multi-function, see page 12)
- 8 Selecting the operating mode
➔ MODE key

9 LEDs for active operating mode

- **PUMP rpm**
Pumping by speed
60 – 6000 rpm
- **PUMP Flow rate**
Pumping by flow rate
- **DISP Time**
Dispensing by time
- **DISP Volume**
Dispensing by volume
- **PAUSE Time**
Pause for intermittent
dispensing
- **PROGRAM**
Program selection 1–4
- **TOTAL**
Read-out of totally delivered
volume

➔ The drive runs only in
clock-wise direction.

⚠ Never run the pump-head
under dry conditions.

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
Control keys

- a RUN/STOP**
Starts and stops the pump
- b MODE**
Changes between operating modes (see page 11)
- c CAL**
Key for saving a set value
- d MAX/RESET**
Multi-function key for:
➔ max. speed (when pump is running)
➔ Reset button resets the calibration to the standard values (default) (see page 23)
➔ Reset button for accumulated volume in Mode TOTAL
- e Increment/Decrement keys**
▲ = increase value
▼ = decrease value
Maintaining pressure on ▲ or ▼ key changes display read-out into fast mode.

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5 Pro4

6

PUMP	*rpm	o Flow rate
DISP	o Time	o Volume
PAUSE	o Time	o PROGRAM
TOTAL	o	

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Start-up information

After switching on the power supply switch, the following values are displayed:

- 1 LED test »8.8.8.8.«
- 2 Firmware version
- 3 Baud rate
- 4 Identification Code (ID-Code) of the pump-head entered in the program last used (e.g.120)
- 5 Program last used
- 6 Displays the operating mode last used

⚠ Before starting the pump the first time or after mounting a new pump-head, the ID-Code must be entered in each program used with this pump-head (see pages 15/44).

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Resetting the parameters of the currently used program to the default values.

- 1 Power switch »OFF«
- 2 Push the **MODE key** whilst switching the drive on

The following parameters are reset:

Mode:	PUMP rpm
Speed:	1000 rpm
Dispensing time:	4.5 sec.
Dispensing volume and flow rate:	depending on pump-head
Pause time:	2.0 sec.

Number of dispensing cycles: 12

➔ The pump-head entered in the currently set program remains stored.

Resetting the parameters of all 4 programs to the default values

- 1 Power switch »OFF«
- 2 Push the **▼ key** whilst switching the drive on.

The following parameters are reset (in all 4 programs):
Pump-head: Model 120
as well as all parameters mentioned above

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Entering the pump-head identification number

In order to obtain correct pumping and dispensing values the individual identification code (ID-Code) of the mounted pump-head must be entered (**⚠ for each of the 4 programs**), **especially** when a new pump is **used for the first time** and each time the pump-head is changed.

- 1 Power switch »OFF«
- 2 Keep CAL-key pressed and
switch power »ON«
- 3 Press ▲ or ▼ key till the
correct ID-Code of the
mounted pump-head is set.
The ID code is indicated on
the pump-head or in the table
on page 44.
- 4 Enter the correct ID-Code with
the CAL-key
(MODE display switches
automatically to PUMP rpm)
- 5 Change to PROGRAM with
the MODE key and select the
next program with the CAL
and ▲▼ keys (confirm with
CAL).


➔ Repeat procedure 1–4 till the ID code of the mounted pump-head is entered in all 4 programs.

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Program selection

When switching the pump on, it always selects the previously used program.

- 1 Change mode to PROGRAM by using the MODE-key.
- 2 Press CAL button, currently set program starts to blink Change to the required program with the **▲▼** buttons
- 3 Confirm with the CAL-key
- 4 The pump automatically returns to the last entered operating parameters of the selected program.

⚠ From now on, any changes carried out in the operating modes are automatically stored in the currently selected program.

TOTAL Volume

Change with MODE key to TOTAL. The totally delivered volume is displayed.

Three display readings are available:

Microlitre: e.g. 17.0 µl = **17.0µ**
Millilitre: e.g. 1700 ml = **17.00**
Litre: e.g. 17.0 l = **17.0L**
(Display in litres from 9999 ml)

For zero-setting, either

- press the reset key twice, or
- switch the pump off.

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
Pumping by drive speed

- 1 Change mode to PUMP rpm (60–6000 rpm) digitally adjustable in stops of 1 rpm
- 2 – Enter required speed with ▲▼ keys
– Start the pump with RUN/STOP
➔ The speed can also be adjusted while the pump is running.
- 3 Maintain pressure on MAX-key while the pump is running
➔ Pump is running at maximum speed (ideal for fast filling or emptying the system)

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Pumping by flow rate

- 1 Change mode to PUMP Flow rate
- 2 Enter the required flow rate with ▲▼ keys (is displayed in ml/min)

Start pump with RUN/STOP

- ➔ The flow rate can also be adjusted while the pump is running.
- ➔ For an accurate flow rate we recommend to calibrate the pump (see page 19).
- ➔ During the pumping process the following values can be retrieved with the MODE-key:
 - Rotation speed (PUMP rpm)
 - Totally pumped volume (TOTAL)

- 3 Maintain pressure on MAX-key while the pump is running
 - ➔ Pump is running at maximum speed (ideal for fast filling or emptying the system)

⚠ Entering the ID-Code of the currently mounted pump-head (pages 15/44) allows the user to set the flow rate in ml/min. This, however, is an approximate value as not yet calibrated.

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



1 Change mode to PUMP Flow rate
Enter the required flow rate with the ▲▼ keys

2 Change mode to DISP Time
– Enter 60 seconds by using the ▲▼ keys
– Start pump with RUN/STOP
➔ Pump stops automatically after 60 seconds
➔ Measure the dispensed liquid by volume or weight and calibrate the ascertained value as follows:

3 Return to mode PUMP Flow rate

4 Press the CAL-button (displayed value blinks)
Enter the weighed or measured value with the ▲▼ keys

5 Confirm with the CAL-key (the flow rate setting returns automatically to the initially entered set point)

6 Start with the RUN/STOP-key

➔ Depending on the application conditions repeated calibration may be necessary.

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
Dispensing by time

The dispensing time can be entered from 0.1s to 999h.

- Change mode to DISP Time
- Enter the required dispensing time with the ▲▼ keys
 ➔ The display accelerates when pressure on the ▲ or ▼ key is maintained.
 Seconds: 0.1–899.9 (in 0.1s steps)
 Minutes: 15'–899' (in 1min steps)
 Hours: 15h–999h (in 1h steps)

- Start pump with RUN/STOP

➔ With the ▲▼ keys the dispensing time can be changed even during the dispensing process. The newly entered set point only takes effect from the subsequent dispensing cycle.

➔ During the dispensing process the following values can be retrieved by pressing the MODE-key:

- speed (PUMP rpm)
- flow rate (PUMP Flow rate)
- totally dispensed volume (TOTAL)

➔ With the ▲▼ keys the speed and flow rate can be changed even during the dispensing process.

➔ Depending on the pump-head and the application, a very short dispensing time can result in dispensing volumes which are not reproducible.

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
Dispensing by volume

- 1 Change mode to DISP Volume
- 2 Use the ▲▼ keys for entering the required dispensing volume. Three display readings are available:
 Microlitre: e.g. 17.0 µl = **17.0** µ
 Millilitre: e.g. 17.00 ml = **17.00**
 Litre: e.g. 17.0 l = **17.0L**
 (Display in litres from 9999 ml)
 ➔ For dispensing accurately we recommend you to calibrate the pump (see p. 22)
- 3 Start pump with RUN/STOP. The pump reduces the rotation speed shortly before the end of the dispensing cycle providing controllable and drop-precise dispensing cycles

➔ The dispensing speed can be adjusted in the modes PUMP rpm and PUMP flow rate.

➔ With the ▲▼ keys the dispensing volume can also be changed even during the dispensing process. The newly entered set point appears shortly in the display taking effect, however, only with the subsequent dispensing step.

➔ During the dispensing process the following values can be retrieved with the MODE-key:

- Rotation speed (PUMP rpm)
- Flow rate (PUMP Flow rate)
- totally dispensed volume (TOTAL)

➔ With the ▲▼ keys the rotation speed and flow rate can be changed even during the dispensing process.

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
Calibrating the volume

- Change MODE to DISP Volume
 ➔ Use the ▲▼ keys for entering the required dispensing volume
- Start with RUN/STOP
 ➔ Pump stops automatically
 ➔ The pump reduces the rotation speed shortly before the end of the dispensing cycle providing controllable and drop-precise dispensing volumes
 ➔ Measure the dispensed liquid by volume or weight
- Stay in the same mode (DISP Volume) and press the CAL-key (display blinks)
- Enter the effectively dispensed volume with the ▲▼ keys
- Save setting with CAL-key
 ➔ Pump is calibrated and returns automatically to the initially entered set point
- Start pump with RUN/STOP
 ➔ Depending on the application conditions repeated calibration may be necessary.

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➔ At the same time, the volume is also reset to the default value

Calibrant reset of volume

- 1 Change mode to DISP Volume
 - 2 Press the CAL button (displayed value blinks)
 - 3 Press the MAX/RESET button
The volume is reset to the default value (Default value blinks).
 - 4 Confirm with the CAL button
- ➔ At the same time, the flow rate is also reset to the default value

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Dispensing by volume within a pre-set time

1 Change to mode DISP Time. Enter the required time with the ▲▼ keys (0.1s – 999h)

2 Change mode to DISP Volume. Enter the required volume in ml with the ▲▼ keys

3 Push the MODE-button
➔ Both LEDs light up simultaneously

4 Start pump with RUN/STOP
➔ The pump dispenses the required volume within the pre-set time.

a **Changes / Corrections**
Press the Reset key
➔ For changing the dispensing time or dispensing volume, MODE returns directly to DISP Time
➔ Repeat steps 1 to 3

b **Error message**
LEDs are blinking alternately

9999 = volume too large
time too short

1111 = volume too small
time too long

Despite volume and time adjustments, the LEDs keep on blinking:
➔ change the pump-head model in order to obtain other flow rate

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Intermittent dispensing (by time)

Intermittent dispensing by time with a pre-set pause time

- 1 Go to mode PAUSE Time
Enter the pause time (between 0.1s and 999h) with ▲▼ keys
- 2 Press MODE key once
➔ Both LEDs for PAUSE and DISP Time light up simultaneously
Enter the dispensing time (0.1s to 999h) with ▲▼ keys
- 3 Start pump with RUN/STOP key
➔ Pump dispenses at intervals until the RUN/STOP button is pushed again.
➔ Active mode blinks.

➔ With the ▲▼ keys the dispensing time can be changed even during the dispensing process. The newly entered set point appears shortly on the display taking effect, however, only with the subsequent dispensing cycle.

- ➔ During the dispensing process the following values can be retrieved with the MODE-key:
- rotation speed (PUMP rpm)*
 - flow rate (PUMP Flow rate)*
 - totally dispensed volume (TOTAL)
- *The rotation speed and flow rate can be changed with the ▲▼ keys.

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1 **MODE**

PUMP	◦ rpm	◦ Flow rate
DISP	◦ Time	◦ Volume
PAUSE	★ Time	◦ PROGRAM
TOTAL	◦	

▲▼

5.0

2 **MODE**

PUMP	◦ rpm	◦ Flow rate
DISP	★ Time	◦ Volume
PAUSE	★ Time	◦ PROGRAM
TOTAL	◦	

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
Please note

For entering the number of dispensing cycles see page 27.

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Intermittent dispensing (by volume)

Intermittent dispensing by volume with a pre-set pause time

- Change mode to PAUSE Time
Enter pause time (between 0.1s and 999h) with the ▲▼ keys
- Press MODE-key twice
➔ Both LEDs for PAUSE and DISP Volume light up simultaneously
Enter the volume with the ▲▼ keys
- Start pump with the RUN/STOP key
➔ Pump dispenses at intervals until the RUN/STOP button is pushed again.
➔ Active mode blinks.

➔ With the ▲▼ keys the dispensing volume can be changed even during the dispensing process. The newly entered set point appears shortly on the display taking effect, however, only with the next dispensing cycle.
➔ During the dispensing process the following values can be retrieved with the MODE-key:
– rotation speed (PUMP rpm)*
– flow rate (PUMP Flow rate)*
– totally dispensed volume (TOTAL)
*The rotation speed and flow rate can be changed with the ▲▼ keys.

If necessary, calibrate the volume as indicated on page 22.

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Number of dispensing cycles

The number of dispensing cycles can be entered when dispensing at intervals (by time and volume).

- 1 Change into one of the two modes for dispensing at intervals
 - DISP Time/Pause, or
 - DISP Volume/Pause
- 2 Push the CAL-key; the number 0 for an infinite number of dispensing cycles blinks
 - ➔ with the RESET-key you can switch between 0 and the number of repetitions last entered
- 3 With the ▲▼ keys the number of dispensing cycles can be adjusted
- 4 Confirm with the CAL-key
- 5 Start dispensing at intervals with the RUN/STOP-key

➔ During the pauses the display shows the remaining number of dispensing cycles.

➔ During long pauses, the remaining interval time can be checked with the CAL-key.

➔ At any time, the dispensing process can be interrupted by pushing the RUN/STOP-key; when resuming the process, the pump continues with the remaining amount of dispensing cycles.

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Dispensing via valve

In order to prevent a **syphoning effect** (liquid flows back to the reservoir) at the moment the pump is stopped, we recommend the use of our MCP-Z *Standard* valve.

This valve can be used for all operating functions which are controllable by the RUN/STOP key.

You can hook up the valve on the right side of the pump and connect it to the corresponding interface (Valves) on the rear panel of the MCP-Z.

When the valve is connected, it automatically opens and closes when the pump is started and stopped.

Start/Stop can also be triggered by an external trigger device, such as

– foot switch (IS 10039)

For further details about the accessories refer to pages 46–48.

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



Dispensing via valve

Mounting the valve

- On the right side of the MCP-Z *Standard* there are two holes for connecting the valve.
- Fix the valve plate by inserting the two screws into the holes. Press plate gently against the drive casing and push it right to the bottom of the hook hole.
⚠ Do not press on the valve.
- Connect the interface cable to the »Valves« socket on the rear panel of the MCP-Z *Standard* pump.
- Insert the tube into the valve as shown on the left.

⚠ If using a new tubing, its hardness, diameter and wall thickness must correspond with the originally supplied tubing:

Silicone tubing
3.2 mm i.d. and 1.6 mm wall,
Order No. MF0037

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
Pumping against pressure

The max. differential pressure depends on the mounted pump-head model (see page 44).

Pump-head with bypass

Some gear pump-head models have an internal fluid bypass.

With the adjusting screw the required differential pressure can be limited to a value between 0.7 bar (10 psi) and the maximum permitted differential pressure of the pump-head used.

Before switching the pump off

⚠ If no non-return valve is mounted on the discharge side of the pump, a siphoning may occur at the moment the pump drive is stopped!
(See page 28–29)

If the pump will not be used for some time, it is strongly recommended to thoroughly rinse the pump-head with water.

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Overcurrent protector

The drive MCP-Z *Standard* features an overload protector. When an overload condition occurs, it is indicated in the display by the letters 'OL' and the pump is stopped.

OL **S** = Speed set-point cannot be reached
OL **I** = Current Overload

Whenever this situation occurs, the pump must be switched off immediately. Before the pump is re-started, it is most important to check the reason for the overload (eg. too high differential pressure, etc.).

➔ Only when the cause of the overload has been detected and the failure corrected accordingly may the pump be started again.

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
Digitale Eingänge (TTL-Pegel)

Digital inputs (TTL-level)

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- Pin 2, local / remote
- Pin 3, stop / run
- Pin 13, speed IN

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Analog interface

Pin 1, GND (ground)

Reference potential for all other inputs

Pin 2, local / remote

For changing between manual control and analog interface. For activating the analog interface, pin 2 must be connected with pin 1 (GND).

Pin 3, stop / run

In remote operation (pin 2 to GND) the pump starts when connected to pin 1 (GND)

Pin 5, speed IN

For external speed control (0–5V, 0–10V, 0–20mA, 4–20mA) Alternatives by means of DIP switches inside the pump (see page 34)

Pin 7, +36V_{DC}

About +36 V_{DC} are available (max. current 1A).

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Analog interface

Pin 9, motor speed

The default setting is $0-10 V_{DC}$, proportionally to the motor speed $60-6000 \text{ rpm}$. Alternatively a frequency range from $0-12 \text{ kHz}$ is available. Alternatives with sliding switch S2 inside the pump (see page 34)

Pin 10, +5 V_{DC}

About $+5 V_{DC}$ are available. (max. current 0.5 A)

Pin 13, speed pin 5 internal

Depending on how the pump is operated, pin 13 has different functions:

Analog interface not activated (Normal operation, i.e. pin 2 is open)
Pin 13 serves as **auto-start function**.

If pin 13 is connected to pin 1 (GND), the pump can be started and stopped directly from the power supply (the power switch must be on).

Analog interface activated (Pin 2 on GND)

- Pin 13 open:
The rotation speed is adjusted via pin 5 (speed IN).
- Pin 13 on GND:
The rotation speed can be adjusted by the speed selector on the control panel of the pump.

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

⚠ Please note

The instrument should only be opened by a qualified technician!

Capacitors inside the pump may still be charged even though the mains plug has been disconnected some time ago.

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Settings of switch S1

Operation via foot-switch

With DIP switch 5 the user can switch between two possibilities:

- »FS toggle« (On/Off)
- »FS direct« (On = as long as pressure is maintained on the foot-switch)

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Pins		DIP-Switch 1	DIP-Switch 2	DIP-Switch 3	DIP-Switch 4	DIP-Switch 5	DIP-Switch 6
Pin 3 (Foot switch)	FS toggle					OFF*	
	FS direct					ON	
Pin 5 speed IN	0–5V		OFF*	OFF*	OFF*		
	0–10V		OFF	OFF	ON		
	0–20mA		OFF	ON	OFF		
	4–20mA		ON	ON	OFF		
Baud Rate RS232	9600 Bd	OFF*					
	1200 Bd	ON					

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Sliding switch S2

This switch affects Pin 9, motor speed

Position A: 0–10 V_{DC} (Standard)
Position B: 0–12 kHz

➔ **Make sure that the pump is disconnected from the mains supply.**

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
Please note

When switching on the MCP-Z *Standard* drive, the serial interface will only respond after a delay of 3 seconds.

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Kommunikationsparameter Communication parameters Paramètres de communication

Standardeinstellung
Default setting
Réglage standard

**9600 baud, 8 bit,
1 stopbit, no parity
(DIP Switch 1 off)**

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Serial interface

RS232 IN (female)
A 9-pin D-socket is available on the rear panel of the pump.

Setting the pump address

Each command string begins with the pump address (1–8). It is factory-set at 1 and can be changed with @. The address will be stored permanently (even after the pump has been switched off). The assignment of the address enables the user to control up to 8 ISMATEC® drives via one interface.

Structure of the commands

The address is followed by a character. Some commands have an additional parameter which always consists of 4 or 5 figures.

The command string **is completed by** a carriage return (ASCII 13). The pump confirms most of the commands with an asterisk *. Yes/No inquiries are answered by + (yes) or – (no).

Multi-digit replies are concluded by (cr) ASCII 13 and (lf) ASCII 10.

Incorrect command strings are answered by #.
Numerical values are confirmed as 3 to 5-digit figures. Four of the five digits are numerals, one digit is either a decimal point or a preceding blank space.

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Each pump requires an individual address
(see page 35).

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Cascading several pumps

Providing that an appropriate software is available, up to 8 ISMATEC® pumps can be controlled via one PC.

In order to establish a cascade you need to connect the pumps with the connecting cable (AG0013) as shown on the left.

The same cable AG0013 serves to connect the first pump to the PC, provided that it has a 9-pin RS232 connector.

If the PC is equipped with a 25-pin connector, you need an additional adapter XC0009.

➔ **Each pump must be allocated its own individual address.**

For this procedure, you first have to connect each pump in turn to the RS232 interface of the PC. Two or more connected pumps must not have the same address. Only when all the pumps have been assigned their addresses can you cascade them as shown on the left.

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Pump-head

In general, the MCP-Z *Process* gear pump is supplied as a complete unit, consisting of the MCP-Z *Process* drive, and as ordered:

- 1 pump-head
- 2 driving magnet
- 3 tubing adapters (2 pcs.)

For the MCP-Z *Process*, a wide range of different gear pump-heads with PTFE, PPS or graphite gears are available.

You find an overview chart of the pump-heads available on page 46.

⚠ Before mounting a new pump-head, make sure that the correct magnet is fitted to the drive.

Mounting the pump-head

Fig. 1

- The following tools are required for mounting a pump-head:
- Allen key $\frac{3}{32}$ "
 - Phillips screwdriver (Size 0)

Remove all screws from the fixing case (a) of the drive.

(Together with the drive, two long screws are additionally supplied for fastening the large pump-heads of the series 200 {at the top and the bottom} of the pump-head face).

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Mounting the pump-head

Fig. 2

Fasten the driving magnet on the motor shaft (through the hole underneath).

➔ Ensure that there is a gap of about 3 mm between the magnet and the brim of the fixing case.

⚠ **During operation the back of the magnet must not touch the housing of the pump drive.**

Fig. 3

Screw the two tubing adapters into the pump-head.

➔ Before screwing the tubing adapters into the pump-head, we recommend you to wrap **PTFE tape around the thread of the adapter**. This ensures optimum sealing and helps to protect the thread.

Fig. 4

Screw the pump-head to the fixing-case of the drive. The face of the pump-head with the engraved model number must be directed to the user. The engraved arrow indicates the flow direction.

➔ Please note that gear pumps can only pump in one direction. An arrow engraved on the front of the pump-head indicates the correct flow direction.

⚠ **Fill the pump-head with liquid before you start pumping. Avoid dry-running of the pump-head.**

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Fig. 2

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For temperatures above 54°C (130°F) consult your ISMATEC distributor for proper seal selection.

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The bypass should only be used for safety reasons and not for pressure control.

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The pump-heads are designed for a system pressure of max. 20 bar (290 psi). The housing is stainless steel 316/DIN 1.4571.

Pump-heads are interchangeable. Please check that the correct magnet is installed.

Please indicate on your order at what temperature you intend to operate the pump.

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

Contain the consummable parts (bushes, seals, gears)

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Accessories

Valve unit

(opens and closes simultaneously with Run/Stop of the pump)

Order No. IS 2047

Non-return valve (Polypropylene)

(included with the drive)

➔ should always be mounted in a vertical position

Order No. SD 0007

Tubing clamps

included with the drive (2 pcs.)

Order No. BN 1001

Tubing adapter 01 for all gear pump-heads

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
Programming Software

Software ProgEdit

- Create your individual application program.
- Download the program via RS232 interface to your MCP Process pump drive.
- Now run your application PC-independently wherever needed.

System requirements:

3 1/2" floppy drive, Win95 or higher, Pentium or higher.

Order No. SOF104

Software Labworldsoft®

is an innovative Windows software for controlling, measuring and evaluating operations in connection with a number of different laboratory instruments.

Hardware requirements

486 PC, at least 33 MHz and 8 MB RAM, mouse, VGA display, monochrome, with at least 16 grey scale values

Software requirement

Windows 3.1 or higher, Windows 95

Order No. SOF103

Software accessories

Extension cable 9F / 9M

9 pins, female / 9 pins, male (2 m long)

Order No. AG0013
(Included with the delivery of the software)

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
Accessories

Interface card Moxa C168

including a 25 pin multiple cable (male) for controlling up to 8 laboratory instruments

Order No. ISA IM0015C168
Order No. PCI IM0030

Adaptor 25 F / 9M

25 pins, female / 9 pins, male for IM0015C168

(1 adaptor for each controlled instrument)
Order No. XC0009

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Maintenance

We strongly recommend you to thoroughly flush the pump-head with water after use.

Provided the MCP-Z *Standard* is operated properly and in compliance with this manual, the moving parts such as gears, etc. are the only parts that are subject to wear and tear.

The Service Kit corresponding to the individual pump-head allows the user to replace the parts that are subject to wear and tear by himself (see page 46).

Repairs

For repairs please send the MCP-Z *Standard* drive to your ISMATEC® representative.

Please give information on defect, date of purchase, serial-no., model.

Spare-parts

After the warranty period your ISMATEC® representative will be pleased to send you on request:

- spare-parts
- parts lists
- wiring diagrams

Please give information on defect, date of purchase, serial-no., model.

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Error check-list

Pump is not running.

Is the pump connected to the mains?
Is the voltage selector set correctly?
Are the correct fuses inserted?
➔ see pages 8–9

Flow rate is incorrect.

Is the correct pump-head ID entered in the actual program?

Have you calibrated the flow rate or dispensing volume in the actual program?
➔ see pages 16/19/22

No response from the RS232 interface?

Is the pump address correct?
Is the pump connected to the correct PC port (COM1,..)?
Do the transmission parameter of the pump and the PC correspond with each other (Baudrate, ...)?
➔ see pages 35–41

Cascaded pumps do not respond.

Are the transmission parameters of the pump and PC correct (Baudrate, ...)?
Have all the connected pumps the correct address?
➔ see page 41

Should problems arise that are not listed here, please contact your local ISMATEC® distributor.

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△ Please note

We also recommend you to observe our **Warranty Terms as well as our Terms and Conditions of Sale**.

In case of any queries, please contact your local ISMATEC® representative.

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Technical Specifications

Drive

Motor type DC-Motor
Speed range 60–6000 rpm
digitally adjustable in steps
of 1 rpm

Differential pressure

max. 6 bar (87 psi) depending on
pump-head see page 44

Remote control

via RS232 and analog interface
(see page 32–41)

Mains connection

230 V_{AC} 50 Hz
115 V_{AC} 60 Hz

Fuse rating

(230 V) 2 x 1.25 A (slow-blow)
(115 V) 2 x 2.5 A (slow-blow)

Protection rating

IP 30

Power consumption

150 Watts max.

Operating conditions

Temperature +5 to +40°C,
Rel. humidity max. 80%
– not condensing, at normal
laboratory conditions

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
Ask for the latest ISMATEC-catalogue.

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ISMATEC® - Your competent partner for demanding metering and dispensing applications.

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A Unit of IDEX Corporation
Laboratoriumstechnik GmbH