



DOWNLOAD
OUR CATALOGUE
www.lamyrheology.com

CONTACT US:
contact@lamyrheology.com
Phone: +33 (0)4 78 08 54 06
Fax: +33 (0)4 78 08 69 44

VISCOMETERS
RHEOMETERS
TEXTURE ANALYZER
TEMPERATURE CONTROL
MEASURING GEOMETRIES



Catalogue 2015-2016

Tactile instruments

PHOTO: BRAINSTORMING.FR Photos: Shutterstock

Dynamic viscosity: η (Eta)

It is defined by the NEWTON equation, and quantify measurement of internal friction of fluid. His determination needs to apply to the fluid a Shear rate (D), and to measure the resistant Shear stress (τ) to this rotation.

Shear rate: D (γ)

is the shearing which subjected by the product in the application. It is known for measurement geometries with small gap. It is not the speed of rotation of the bob (in rpm !).

Either a sheared fluid, by a laminar move (dV), between two parallel plates with a surface (S) and separate by a distance dx.

Shear stress: τ (Tau)

There is the shearing force (F), with which the sample answers to the shear rate (D), divided by the contact surface (S).

Rheology:

There is the « science » of « flow ».

Associated physical measurements, realised with the hand of Rheometers, enables the visualisation of the behaviour of the product in various flow, temperature and time conditions.

Rotating viscometer:


a - With coaxial cylinders

The fluid is sheared between two coaxial cylinders, with radius Ri and Ra and a length L, by a laminar move which are breaking down in multi-layer with different angular speed from 0 (for the layer in contact with the fixed cylinder) to ω (for the layer in contact with the rotating bob). The relative move of layers towards others give, a shear rate D and one Shear stress τ .

By imposing ω and measuring M, the resisting torque to this rotation, we calculate D and τ according :

$$\delta = Ra / Ri \quad Ri / Ra \rightarrow 0.92$$

Shear stress:
 $T_{rep} = (1 + \delta^2 / 2 \delta^2) * (M / 2\pi L R_i^2)$
 Shear rate:
 $D_{rep} = \omega * (1 + \delta 2) / (\delta 2 - 1)$




Rq : The determination of D is possible only if the gap is small. (i.e. DIN / ISO 3219 Standard).

$$\tau = \eta * D \text{ in Pa.s}$$

For memory:

$$1 \text{ Pa.s} = 10 \text{ Poises or } 1 \text{ mPa.s} = 1 \text{ cPoises}$$



$$D = dV / dx \text{ in s}^{-1}$$

$$\tau = F / S \text{ in Pa (N / m}^2\text{)}$$


Rheograms:

displayed curves of the flow behaviour of a fluid.

The curves $\tau = f(D)$ enables, by adapted fitting, the access to direct related parameters with the application.

b- With Cone-Plate :

The fluid is placed between a Plate and a Cone with angle δ ($< 3^\circ$). The cone, maintained to a constant speed induce a laminar shearing move. In those conditions, τ and D are constant in the gap, according :



Shear stress / Shear rate

$$\tau = 3M / 2\pi R^2 / D = \omega / \text{arc } \delta$$

Rq : You would be vigilant on the sample volume including in the gap, because the great influence of the radius R on the τ value !


Study of different rheological behaviours

TYPE	NEWTONIAN	PSEUDO-PLASTIC	PLASTIC
Description	A sample is named newtonian when his viscosity stay constant, whatever the shear rate is. It's not necessary to define exactly the shear rate for the measurement, just the temperature is important.	One pseudo-plastic sample has a viscosity which decrease when the shear rate increase. This flow behaviour is due to the molecules form and to their internal structure.	One sample presents a plastic behaviour, when his viscosity decrease when the shear rate increase, but from one original shear stress upper than 0, called YIELD VALUE (τ_{0D}), shear stress under which the product doesn't flow. It behave like a solid body.
Rheogram			
Viscosity			
Equations	$\eta = \tau / \alpha$	$\tau = K * D^n$ ou $n < 1$	$\tau = \tau_0 + \eta_{lo} * D$
Examples	<ul style="list-style-type: none"> Water: 1 mPa.s to 20° C Oils: 150 to 400 mPa.s (motor) 300 to 800 mPa.s (gears) Mercury: 1.5 mPa.s Gas: 0.01 to 0.02 mPa.s 	<ul style="list-style-type: none"> Coating, Varnish, Cosmetics, Mineral Suspensions... 	<ul style="list-style-type: none"> Toothpaste, Ointment, Grease, All very concentrated suspensions...

The thixotropy

One thixotropic product is a sample for which the variation of viscosity in function of shear rate is associated to a variation through the time.

Owe talk about Thixotropy or Rheopexy, with the condition of REVERSIBLE Transformations: frozen or solidification.



Causes of thixotropy :

- Molecular structure
- « Château de cartes » with layers
- Particules mixing
- Ball loose Package...

CUSTOMER SERVICES

GLOSSARY

Adhesiveness: is the sticky power of a product. It is measured during a tension phase in texture analysis, by the negative force measured and also by the surface under the base line.

ASTM: American Society of Testing Materials. American organisation in charge of creating ASTM standards.

BINGHAM: model of rheological flow behaviour, characteristic of plastic products (shear-thinning with yield stress).

CASSON: model of rheological flow behaviour, allows the precise determination of non-linear plastic product's yield stress.

Centipoise (cP): measuring unit of dynamic viscosity in the MKSA system; equivalent to mPa·s in the SI system.

Coaxial cylinders: one cylinder with cap contains the product (cup) and one cylinder of a smaller size and another cylinder rotates inside (measuring bob) and imposes shear rate (D) known in the sample. (see DIN Standard).

Cone-plate: measuring geometry composed of one plate on which the product to be measured is placed and a low-angle cone (2° max), which shears the sample.

Consistency: notion of force with which a product resists compression. Quantified in texture analysis by Maximum Force that is measured during a compression phase.

Coquette principle: principle of rheometer function in which the cup or the lower plate turns or oscillates, and the measuring bob or cut or upper plate measures torque. This principle lets you separate the part deforming the sample from the part that measures.

D (or $\dot{\gamma}$): shear rate actually subject to the fluid to be measured, expressed as s⁻¹.

Dilatancy: increase of viscosity with the effect of rotation speed.

DIN: German Original Standard, specifying measuring geometries at a defined shear rate. Became ISO 3219.

Elasticity: Ability of a sample to recover its initial state after having been deformed. Inversely proportional to the relaxation % in texture analysis tests.

ETA (η): Dynamic Viscosity: quantifies a fluid's internal frictions; determined by the rotating principle; torque resistant to rotation; expressed in Pa·s.

K: consistency coefficient according to the Ostwald model; it shows a product's viscosity at 1 s⁻¹.

KREBBS Unit: viscosity measuring unit obtained with a KU110 measuring bob, at 200 rpm.

M (mNm): measured torque in response to the rotation of the measuring bob based on the product's viscosity.

Measuring bob (spindle): element immersed in fluid which rotates and measures the resistant torque of a product, according to the Searle principle.

Measuring geometry: set of spindles and cups or cones and plates used to measure viscosity. It enables, if well defined, to control the shear rate (D) subjected by the product.

N: rotation speed of motor, in rpm, which generates a shear rate (D) which depends on the measuring geometry used.

n: behaviour index of the Ostwald model; shows shear-thinning character of a product.

NEWTON: model of rheological behaviour model, characterising fluids for which only temperature has an influence on viscosity.



DISTRIBUTORS



EUROPE

Belgium - Bosnia - Bulgaria
Croatia - France - Germany
Greece - Hungary - Ireland
Italy - Netherlands - Poland
Romania - Serbia - Spain
Ukraine - United Kingdom

AMERICAS

Brazil - Canada
Colombia - Guatemala
Honduras - Mexico
Salvador
United States

AFRICA

Benin - Egypt
Nigeria - Tunisia

ASIA

China - India - Indonesia
Iran - Israel - Japan - Jordan
Lebanon - Malaysia - Pakistan
Philippines - Qatar - Russia
Singapore - South Korea
Taiwan - Turkey
United Arab Emirates
Vietnam

OCEANIA

Australia
New Zealand

Go to www.lamytheology.com
to see a list of our distributors