


Degassing a microfluidic system

Dissolved gas and trapped air bubbles in a microfluidic system can affect repeatability and stability. Degassing a fluid system is easy using the LabSmith Breadboard Reservoir (BBRES-C360 or BBRES-T116) or Chip Reservoir (C360-405R)


Using a breadboard reservoir:

- construct your fluid circuit as desired using the breadboard reservoir. Make sure the circuit is closed (plug any open ports)
- fill the reservoir at least half-full with the desired fluid
- secure the cap on the reservoir with an o-ring seal
- connect a manual luer-lock syringe to the reservoir cap with a short piece of tubing
- use the syringe to pull a vacuum in the fluid circuit



Using a chip reservoir:


- connect a chip reservoir to a bonded port connector on a microfluidic chip
- ensure all open ports on the chip are plugged
- fill the reservoir at least half-full with the desired fluid
- connect a manual luer-tip syringe to the chip reservoir and pull a vacuum



Degassing a microfluidic system using the LabSmith Breadboard Reservoir

Dissolved gas and trapped air bubbles in a microfluidic system can affect experiment repeatability and stability.

- construct your fluid circuit as desired using the breadboard reservoir (BBRES-C360 or BBRES-T116).
- ensure the circuit is closed (plug any open ports)
- fill the reservoir at least half-full with the desired fluid
- secure the cap on the reservoir with an O-ring seal
- connect a manual luer-lock syringe to the reservoir cap with a short piece of tubing
- use the syringe to pull a vacuum in the fluid circuit



Degassing Tips

- Degas the circuit with syringes in the dispensed position to make it easier to remove the air bubbles
- Put a small amount of fluid in the manual syringe so you can observe when the fluid degassing is complete (when you no longer see bubbles rising into the fluid)
- If you continue to observe bubbles in the fluid you may have a leak in your fluid circuit or a blocked path (see Degassing Tips page)
- A ~10 ml syringe is a good size for degassing a circuit, however several iterations may be necessary to removed all trapped gas.
- Sonicate your fluid prior to introducing it into your system or reservoir to remove most of the dissolved gasses

Degassing Troubleshooting

If you continue to observe air bubbles in your system after you have degassed you may have a leaky or blocked connection

Problem	Solution
Incorrectly assembled one-piece fitting	Tubing must extend past the end of the one-piece fitting to properly seal
Incorrectly cut tubing	Clean straight cut required for leak-free installation
Pinched tubing	Cut PEEK tubing with a sharp blade to minimize risk of pinching ends, look for kinks along the tubing
Blocked port or valve	<ul style="list-style-type: none">• Always flush components (especially valves) after use if working with particles or condensing solutions• If using fused-silica tubing, incorrectly cut ends can cause glass particles to break off and block ports
Broken o-ring on reservoir cap	Check to ensure the o-ring isn't pinched or broken