

BOD analyzer with Biosensor, Model Quick BOD α1000

Specifications	
Model	Quick BOD α 1000
Measurement method	Biosensor
DO electrode	Clark cell type oxygen electrode
Measurement range	2 ~ 50 mg/L (basic measurement range based on trichosporon cutaneum membrane)
Measurement time	Approx. 60 minutes per measurement
Measurable sample	Liquid samples containing soluble organic
Calibration	1 point calibration by glucose glutamic solution
Repeatability	±5% (in case of 20-50 mg/L F.S. setting, using trichosporon cutaneum membrane and stable environment temperature)
Temp. control	20 to 40 degree C (adjustable) Automatic warming/cooling PI control function with peltier device
Output	Print-out output BOD value hold output (DC 0-5 V) Probe output (DC 0-5 V)
Power requirement	100VAC, 50/60 Hz, 3 A
Dimensions	260(W) × 320(D) × 409(H) mm
Weight	Approx. 16 kg, not including reagents
Optional accessories	Microbe membrane (trichosporon cutaneum) Pipettes or measuring apparatus



Order No.	Descriptions	Remarks	Quantity
C00018601	BOD Analyzer Quick BOD α 1000		
Standard Accessories	Printer	● Printy 3	1
	Microbial immobilization kit	● Immobilization membrane kit (10 sets/pkg) ● Immobilization tool (1 pair)	1
	BOD seeds	● 5 capsules	1
	DO probe kit	● DO electrode (1) ● Membranes (5) ● Electrolyte (1) ● Metal polish (1) ● Syringe (1)	1
	0.5M Phosphate buffer solution	● 5L	1
	5000 mg/L BOD standard solution	● 500mL	1
	Containers	● Plastic vial, 30 mL (1) ● Plastic bottle, PP, narrow-neck, 500 mL (1) ● Plastic bottle, PP, wide-neck, 100 mL (2) ● Plastic bottle, PP, wide-neck, 500 mL (2) ● Poly flattened container, wide-mouth, 10 L (2) ● Poly flattened container, wide-mouth, 20 L (1)	1

Optional accessories			
C00041303	Microbe membranes, 5/pkg	● Trichosporon cutaneum	1

[Apparatus for the estimation of biochemical oxygen demand (BODs) with microbial sensor.]

JIS K 3602-1990 (established Sept. 1st, 1990)

Although the BODs (JIS K 3602) microbe electrode method and BOD5 (JIS K 0102) 5-day method are different, the expected BOD5 value can be approximated from the BODs value. This enables water quality to be monitored over a short period of time in order to permit any necessary measures to be enacted quickly and effectively. It was therefore incorporated as a separate method.

☆ This analyzer complies with the JIS standard only if the optional microbe membrane (Trichosporon cutaneum)

※ Subject to change

ISO 9001 certified



CENTRAL KAGAKU CORP.

Headquarters:
Shouei Bldg., Hongo 3-23-14,
Bunkyo-ku, Tokyo 113-0033, Japan

Telephone : +81-3-3812-9186
Facsimile : +81-3-3814-7538

Tokyo, Osaka, Nagoya, Fukuoka

Distributed by

BOD analyzer with Biosensor

Quick BOD α1000

- Highly correlated with the BOD₅ method by microbe immobilization
- Quick measurement time (60 min.)




Quick BOD α 1000 is a BOD analyzer with biosensor for BOD measurement equipping a high precise liquid transition, temperature control and detection. Microbe immobilization kit enables to measure BOD values by using microbes in customer's activated sludge optimized for customer's sample in addition to normal microbes. Furthermore the detection by Clark cell type DO probe provides long stable measurements.




Highly correlated with the BOD₅ method

By using the membrane on which microbe was immobilized with microbe immobilization kit as standard accessory, you can get a high correlation with the BOD₅ method. See the figures below.



[Procedure of preparing Immobilized microbe membrane]



The biosensor method is the fastest!!

Normally BOD measurements require 5 days. The analyzer can measure BOD in 60 minutes due to Biosensor technique and the incubation under 20 degree C temperature for 5 days periods according to the standard method, No. JIS K 0102, is not required.

Touch key pads easy to operate, compact size for the minimum space.

The analyzer features a simple menu structure and easy touch keys to operate. The analyzer has a slim and compact body which stands upright and does not require large space to install.



Easily viewable display (20 digits x 4 lines, back-light LCD). Simple operating keys.


**Basic measurement range 2-50 mg/L
Up to 100 mg/L available**

The basic measurement range is 2 mg/L to 50 mg/L. The maximum measurement range is expandable up to 100 mg/L. It is adaptable for monitoring or controlling the water quality of not only effluent but also influent or treated sewage in the process.


[Principle of measurement]

The change in respiratory action caused by the consuming of organic matter by bacteria is monitored by an oxygen electrode. This bacteria is incorporated into the oxygen electrode. When a solution containing organic matter is fed to the microbe electrode sensor, the microbes ingest the organic matter increasing respiratory action and then consuming oxygen. The amount of oxygen passing through the microbe membrane decreases, resulting in a change in the current which is output by dissolved oxygen electrode. This change corresponds to the concentration of organic matter.

[Sample water not containing organic matter]



[Sample water containing organic matter]




Continuous measurement available

The high-precision auto temperature control function equipped newly supports the stable operation. The continuous measurement for waster water monitoring is available.

Low and easy maintenance


The newly equipped components such as Clark cell type oxygen electrode, antibacterial tubing and check valve provide stable measurement and low maintenance.

[Inside layout]



[Comparison between the biosensor method and the BOD₅ method]

**Example 1
Effluent from a chemical plant**



**Example 2
Wastewater from an electric components manufacturer**

