

WE NEED THIS TO BE
Easy to use

NANOSPRAY
IN 5 MINUTES

...BUT WE DON'T WANT TO GIVE UP ANY
Performance



NATHAN YATES
UNIVERSITY OF PITTSBURGH

You asked for it. We listened.

- Plug & spray operation
- Single connection
- Temperature-control
- Factory-verified voltage connections
- Chip-to-chip consistency
- Inexpensive to operate


PICOCHIP.
NANOSPRAY COLUMNS




Order PicoChip in a variety of colors, whether to color-code experiments or to designate users

Model PCH-550
For Thermo IDQ™, Orbitrap™,
Velos™, QExactive™, and TSQ™

The proven high-performance of the PicoFrit® column is integrated with a factory-configured high-voltage liquid-junction connection in a compact, optimized, rugged format. The PicoChip is not only easy to use, but easy on the budget, too. And it's all backed by New Objective's commitment to service and quality.



Place PicoChip on source




Tip is protected inside chip




Insert chip key to expose tip



Tip is exposed to spray



Collect data

 NEW OBJECTIVE

New Objective, Inc.
Two Constitution Way
Woburn, MA 01801-1023

888 220 2998 Orders
781 933 9560 Phone
781 933 9564 Fax

www.newobjective.com
sales@newobjective.com

Innovation in High-Sensitivity LC-MS Research



* Introductory Offer - Cannot be combined with other offers or discounts.
Call New Objective for more information.

PicoChip and PicoFrit are trademarks and registered trademarks of New Objective, Inc. All rights reserved. LTO, Orbitrap, Velos, QExactive, and TSQ are trademarks of Thermo Scientific. PicoChip is manufactured and sold under license of U.S. Patents 6,396,183, 6,190,559, and 5,997,740.