

SICCO Vacuum-Desiccator

Panels:	Material:	Temperature resistance:	Use:	Stand:
crystal clear	PMMA	-20 °C to +70 °C	Vacuum up to 10 ⁻⁴ Torr	rubber feet

Specification:

solid walls of 25 mm thickness, circular silicone gasket, vacuum gauge, V 1880-04 including two shelves, usable with up to three shelves, V 1880-08 and V 1880-12 including each four shelves, usable with up to each eight shelves

A	Cat. No. Vacuum 1:	V 1880-04
	Overall dimensions* W x H x D mm	300 x 350 x 300
	Excess length of valves + vacuum gauge mm	100
	Usable interior space W x H x D mm	250 x 260 x 250
	Usable volume liters, approx.	16
	Total volume liters, approx.	16
	Weight kg/approx.	15
	Maximum all-over load per shelf kg	9
	Total all-over load kg	30
B	Cat. No. Vacuum 2:	V 1880-08
	Overall dimensions* W x H x D mm	300 x 450 x 300
	Excess length of valves + vacuum gauge mm	100
	Usable interior space W x H x D mm	250 x 360 x 250
	Usable volume liters, approx.	22
	Total volume liters, approx.	22
	Weight kg/approx.	18
	Maximum all-over load per shelf kg	9
	Total all-over load kg	40
C	Cat. No. Vacuum 3:	V 1880-12
	Overall dimensions* W x H x D mm	310 x 560 x 310
	Excess length of valves + vacuum gauge mm	100
	Usable interior space W x H x D mm	250 x 460 x 250
	Usable volume liters, approx.	29
	Total volume liters, approx.	29
	Weight kg/approx.	27
	Maximum all-over load per shelf kg	9
	Total all-over load kg	50



Add. accessories

Pages 31–34.

*Overall dimensions are indicated without hinges and valves; height is measured starting from base