+1 800 426 0191 +1 360 679 2528 www.idex-hs.com

Valve Wrenches

- ▶ For convenient wrench-tightening of fittings on high pressure rotary shear valves
- ▶ For removal of knobs on Upchurch Scientific[®] Manual Valves

The smartly designed IDEX Wrench is a double-ended slotted socket wrench that fits over 1/16" and 1/8" OD tubing. It easily loosens and tightens 1/4" and 5/16" hex head stainless steel or PEEK fittings. The "Z" shape of the IDEX Wrench provides ideal leverage for changing sample loops and fittings, and keeps one end from restricting the use of the other.


The V-103 is an Allen (hex-key) wrench designed to remove the knob from Upchurch Scientific V-101 valves (page 133). The V-104 is an Allen wrench that can be used to remove the knob from Upchurch Scientific Medium Pressure Selection and Injection Valves (also found on page 133).

MXX Replacement Fittings

Use these replacement Ferrules and O-rings for 1/8" and 1/16" tubing with the MXX Series II valves shown on page 133. Please see the part number chart below for a list of individual part numbers.

Biocompatible

6810 0.050" Allen Wrench (replacement) V-103 V-104 5/64" Allen Wrench (replacement) IG BRACKET ACCESSORIE MOUNT 7160 Mounting Panel Valve Angle Bracket 7160-010 7160-029 Ring Stand Mounting Bracket VALVE BRACKE M-615-1 Mounting Bracket for Upchurch Scientific® Switching Valves M-615-2 Mounting Bracket for Upchurch Scientific Injection and Selection Valves Part No. Description Qty. 7770-039 Ferrules for 1/8" OD Tubing 25-pk 7770-040 Ferrules for 1/8" Tubing 50-pk Ferrules for 1/8" Tubing 100-pk 7770-041 Ferrules for 1/16" OD Tubing 7770-044 25-pk Ferrules for 1/16" Tubing 7770-045 50-pk 7770-046 Ferrules for 1/16" Tubing 100-pk





Rheodyne® mounting brackets and panels of different shapes and

sizes organize and provide a sturdy support for Rheodyne valves.

The Ring Stand Mounting Bracket now allows the valves to mount onto

Mounting Brackets

common laboratory equipment.

