




# SafetyCaps Thread GL 45 with shut-off



## Simplified change of bottles

### No air entrapment in capillaries



After the supply container is changed, the current analysis can be continued without interruption.

### Ready for use again - quickly

During repairs on HPLC pumps, the shut-off on the SafetyCaps can be closed so the capillaries and solvent filter remain in the solvent. Thus the pump can be quickly and easily flushed when returning to operation after repairs.



# GL 45 with shut-off



| Fig.     | Part No. | Description                                   | Thread size  | Connector<br>Ø 3.2 mm OD | of which<br>with shut-off |
|----------|----------|---|--------------|--------------------------|---------------------------|
| <b>A</b> | 107 119  | SafetyCap I with shut-off                     | GL 45        | 1x                       | 1x                        |
| <b>B</b> | 107 919  | SafetyCap II with shut-off                    | GL 45        | 2x                       | 2x                        |
| <b>C</b> | 107 920  | SafetyCap III with shut-off                   | GL 45        | 3x                       | 3x                        |
| <b>D</b> | 107 519  | SafetyCap VI with shut-off                    | GL 45        | 6x                       | 6x                        |
| <b>E</b> | 107 219  | SafetyCap II with shut-off <b>(combined)</b>  | GL 45        | <b>2x</b>                | <b>1x</b>                 |
| <b>F</b> | 107 319  | SafetyCap III with shut-off <b>(combined)</b> | GL 45        | <b>3x</b>                | <b>2x</b>                 |
| <b>G</b> | 117 010  | Air valve for SafetyCaps                      | UNF 1/4" 28G | -                        | -                         |
| <b>H</b> | 197 010  | Air valve for SafetyCaps (10 pc./unit)        | UNF 1/4" 28G | -                        | -                         |