

# SafetyCaps Special threads

**A** 107 005



1x



GL 28

**B** 107 006



2x



GL 28

**C** 107 100



1x



S 40 / GL 40

**D** 107 101



2x



S 40 / GL 40

**E** 107 105



1x



1x





S 40 / GL 40

## Thread S 40




Now you can use our proven safety system without an adapter, directly on supply bottles with thread sizes S 40 and GL 40. Many solvents are already delivered in such containers. With the new sealing system, you can connect the containers directly to the HPLC system, without decanting or thread adapters.





### Include the suitable container to your order

The chapter "containers", starting from page 87, shows all suitable canisters and lab bottles for SafetyCaps.

Fig.	Part No.	Description	 Thread size	 Connector Ø 3.2 mm OD	 of which with shut-off
<b>A</b>	107 005	SafetyCap I (GL 28)	GL 28	1x	-
<b>B</b>	107 006	SafetyCap II (GL 28)	GL 28	2x	-
-	107 512	SafetyCap III (GL 38)	GL 38	3x	-
<b>C</b>	107 100	SafetyCap I (S 40)	S 40 / GL 40	1x	-
<b>D</b>	107 101	SafetyCap II (S 40)	S 40 / GL 40	2x	-
-	107 742	SafetyCap III (S 40)	S 40 / GL 40	3x	-
<b>E</b>	107 105	SafetyCap I (S 40) with shut-off	S 40 / GL 40	1x	1x
-	107 030	SafetyCap I (GLS 80)	GLS 80	1x	-
-	107 031	SafetyCap II (GLS 80)	GLS 80	2x	-
<b>F</b>	107 032	SafetyCap III (GLS 80)	GLS 80	3x	-
<b>G</b>	107 035	SafetyCap IV (B 83)	B 83	4x	-