


SafetyWasteCaps Thread S 55

A 107 917

-  3x
-  •



B 107 924

-  2x
-  1x
-  •



C 107 936

-  3x
-  1x
-  •



D 107 960

-  2x
-  •
-  •



E 108 200

-  2x
-  •
-  •




Helpful ordering tool

An overview, showing all possible connections to the tube connector at SafetyWasteCaps, you will find on the helpful fold out pages at the end of this chapter.

Thread information

The number "51" is often in the cover.

ca. 53.5 mm



What kind of thread my canister has?

An overview, showing all helpful thread informations and measurements, you will find on the pages 87 - 88.













What is included in the scope of delivery?

The shown pictures accord to the scope of delivery. Fittings for capillaries with 2.3 / 3.2 mm outer diameter and if applicable, a tube fitting for tubes with 6.4 - 9.0 mm inner diameter are included.

Protect from harmful vapors!

Due to different container sizes and space conditions on smaller SafetyWasteCaps we offer exhaust filters in various sizes. Just choose your suitable size **S, M and L** for service lifes from 3 up to 6 months. All exhaust filters you will find on the pages 59 - 60.

Fig.	Part No.	 Capillary connector 2.3 / 3.2 mm OD	 Tube connector 6.4 - 9.0 mm ID	 Connector for exhaust filter	 Level control (mechanical)	 Level control (electronic)	 Safety funnel with automatic closure	 Safety funnel with shut-off	 Ground connection
A	107 917	3x	-	●	-	-	-	-	-
B	107 924	2x	1x	●	-	-	-	-	-
C	107 936	3x	-	●	-	-	-	-	●
-	107 943	2x	1x	●	-	-	-	-	●
D	107 960	3x	-	●	●	-	-	-	-
-	107 963	2x	1x	●	●	-	-	-	-
-	108 030	3x	-	●	●	-	-	-	●
E	108 200	3x	-	●	-	●	-	-	-
-	108 201	2x	1x	●	-	●	-	-	-
F	108 033	2x	-	●	-	-	●	-	-
-	108 128	-	-	●	-	-	●	-	-
G	118 143	2x	-	●	-	-	-	●	-
-	118 149	-	-	●	-	-	-	●	-