

Accessories for Autoclave

Printer^{※1}

Accessories: Thermal paper, paper holder and guide
 ※Factory Option



Sterilizing Container for Wastes ST-ZERO

Accommodating 25 × 25ml disposable pipettes
 Accessories: Two stands



Deodorant for Autoclave

300g × 6



Sterilizing Bag

DAB-700

Allowable Temperature Max 180°C : For SX-700 · SX-700E (100 pieces)

DAB-240 (fig.)

Allowable Temperature Max 121°C : For ES-215 (200 pieces)

DAB-321 (fig.)

Allowable Temperature Max 180°C : For SX-500 · SX-500E · BSX-500, SX-300 · SX-300E · ES-315 (100 pieces)

DAB-120

Allowable Temperature Max 140°C : Opaque 4 layered structure supplied in aluminum foil.
 For SR-240 (50pieces)




■ table of Basket, Bucket, and Stand for Testing Durham's Tubes combination


Autoclave	Basket, Bucket, Stand for testing tubes	Model	Inside dimension (diameter × height)	Remarks
ES-215	Stainless basket	ACA-215B	φ 224 × 376mm	(1 pc installable)
	Adjustable shelf	ACA-215S	—	Stainless Steel
	Stainless bucket	ACA-215K	φ 222 × 364mm	(1 pc installable)
	Stand for testing Durham's tubes	STD-215	φ 219 (outside diameter) × 180mm	hole diameter φ 19 × 72 (2pcs. installable)
SX-300 · SX-300E · SX-500 · SX-500E · BSX-500 · ES-315	Stainless basket	ACA-315B	φ 300 × 182mm	(3pcs. installable, 2pcs. installable for SX300/SX-300E)
	Stainless bucket	ACA-315K	φ 299 × 279mm	with look for stacking (2pcs. installable, except SX300/SX-300E)
	Stainless bucket	ACA-500K	φ 299 × 519mm	(except SX300/SX-300E)
	Stand for testing Durham's tubes	STD-500	φ 294 (outside diameter) × 180mm	hole diameter φ 19 × 128 (3pcs. installable, 2pcs. installable for SX300/SX-300E)
SX-300 · SX-300E · SX-500 · SX-500E · BSX-500	Stainless bucket	ACA-300K	φ 299 × 279mm	without hook for stacking
SX-700 · SX-700E	Stainless basket	ACA-700B	φ 345 × 181 mm	(3pcs. installable)
	Stainless bucket	ACA-700K	φ 344 × 279mm	with look for stacking (2pcs. installable)
	Stand for testing Durham's tubes	STD-700	φ 344 (outside diameter) × 180mm	hole diameter φ 19 × 172 (3pcs. installable)
SDL-320	Stainless basket	ACA-315B	φ 300 × 182mm	(3pcs. installable)
	Stainless bucket	ACA-315K	φ 299 × 279mm	(2pcs. installable)
ES-215 · SR-240	Stainless basket	ACA-240B	φ 224 × 181 mm	(2pcs. installable)

■ table of Stainless Long Basket and Lid combination (For Sterilizing Bag)

ES-215	Stainless Long Basket	ACA-215T	φ 224 × 375mm	(Sterilizing Bag DAB-240)
	Lid	215TL	—	(Sterilizing Bag DAB-240)
SX-300 · SX-300E	Stainless Long Basket	ACA-300T	φ 300 × 376mm	(Sterilizing Bag DAB-321)
	Lid	300TL	—	(Sterilizing Bag DAB-321)
SX-500 · SX-500E · BSX-500 · ES-315	Stainless Long Basket	ACA-315T	φ 300 × 557mm	(Sterilizing Bag DAB-321)
	Lid	300TL	—	(Sterilizing Bag DAB-321)
SX-700 · SX-700E	Stainless Long Basket	ACA-700T	φ 345 × 561mm	(Sterilizing Bag DAB-700)
	Lid	700TL	—	(Sterilizing Bag DAB-700)



■ dimension of chamber, Basket, Bucket, and Stand combination



● The image shown for the product may not display its real color or shape.